

Anti-TLR7 antibody (175-225) (STJ503281)

STJ503281

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TLR7 (175-225) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Cat, Human, Mouse, Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.72 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51284
Gene Symbol	TLR7
Uniprot ID	TLR7_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 175-225 on human TLR7 protein.
Immunogen Region	175-225
Specificity	
Immunogen Sequence	



Western blot with Anti-TLR7 Antibody (STJ503281) on HepG2 cells. Antibody at 1:500 dilution in buffer.