

Anti-TLR12 antibody (440-470) {Biotin} (STJ503292)

STJ503292

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-TLR12 (440-470) is suitable for use in ELISA and Western Blot research applications. |
| Applications | ELISA, WB |
| Host/Source | Rabbit |
| Reactivity | Chimpanzee, Human, Monkey |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.55-0.69 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 384059 |
| Gene Symbol | Tlr12 |
| Uniprot ID | TLR12_MOUSE |
| Immunogen | Synthetic peptide taken within amino acid region 440-470 from mouse TLR12 protein. |
| Immunogen Region | 440-470 |
| Specificity | |
| Immunogen Sequence | |