

Anti-TLR13 antibody (400-480) (STJ503294)

STJ503294

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TLR13 (400-480) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	279572
Gene Symbol	Tlr13
Uniprot ID	TLR13_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 400-480 of mouse TLR13 protein.
Immunogen Region	400-480
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat tissue with Anti-TLR13 primary antibody (STJ503294) at 1:100 dilution in buffer. Section was treated with DAB and Haematoxylin stain.