

Anti-TLR5 antibody (195-225) {FITC} (STJ503302)

STJ503302

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TLR5 (195-225) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7100
Gene Symbol	TLR5
Uniprot ID	TLR5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 195-225 on human TLR5 protein.
Immunogen Region	195-225
Specificity	
Immunogen Sequence	