

Anti-TLR9 antibody (430-460) {FITC} (STJ503308)

STJ503308

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TLR9 (430-460) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Chimpanzee, Human, Monkey, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.57-0.76 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	54106
Gene Symbol	TLR9
Uniprot ID	TLR9_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 430-460 on human TLR9 protein.
Immunogen Region	430-460
Specificity	
Immunogen Sequence	