

## Anti-TM6SF2 antibody (325-375) {Biotin} (STJ503316)

STJ503316

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-TM6SF2 (325-375) is suitable for use in Dot Blot, ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	DB, ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.68-0.70 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">107770</a>
<b>Gene Symbol</b>	<a href="#">Tm6sf2</a>
<b>Uniprot ID</b>	<a href="#">TM6S2_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 325-375 on human TM6SF2 protein.
<b>Immunogen Region</b>	325-375
<b>Specificity</b>	
<b>Immunogen Sequence</b>	