

Anti-TMC1 antibody (710-760) {FITC} (STJ503320)

STJ503320

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TMC1 (710-760) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	117531
Gene Symbol	TMC1
Uniprot ID	TMC1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 710-760 on human Transmembrane channel-like protein 1.
Immunogen Region	710-760
Specificity	
Immunogen Sequence	