

Anti-Timeless antibody (1178-1197) {FITC} (STJ503335)

STJ503335

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Timeless (1178-1197) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.98-1.26 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	21853
Gene Symbol	Timeless
Uniprot ID	TIM_MOUSE
Immunogen	Synthetic peptide: SGTPRVHRKKRFQIEDED corresponding to amino acids 1178-1197 on Mouse Timeless protein.
Immunogen Region	1178-1197
Specificity	
Immunogen Sequence	