

Anti-TOMM40 antibody (300-350) {FITC} (STJ503347)

STJ503347

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TOMM40 (300-350) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10452
Gene Symbol	TOMM40
Uniprot ID	TOM40_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on human Mitochondrial import receptor subunit TOM40 homolog.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	