

Anti-TPCN1 antibody (STJ503363)

STJ503363

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TPCN1 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10,000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53373
Gene Symbol	TPCN1
Uniprot ID	TPCN1_HUMAN
Immunogen	Synthetic peptide taken within TMD7 and TMD8 on human TPC1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-TPCN1 Antibody (STJ503363) with TPC1 peptide (A), MCF7 cells (B) and SKNBE2 cells (C). Antibody at 1:500 dilution in buffer.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081