

Anti-TPCN2 antibody (205-255) (STJ503367)

STJ503367

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TPCN2 (205-255) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.57 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	219931
Gene Symbol	TPCN2
Uniprot ID	TPC2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 205-255 on human TPC2 protein.
Immunogen Region	205-255
Specificity	
Immunogen Sequence	