

Anti-TPD52 antibody (1-50) {FITC} (STJ503376)

STJ503376

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TPD52 (1-50) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.69 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:1, 000 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7163
Gene Symbol	TPD52
Uniprot ID	TPD52_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human Prostate leucine zipper isoform 5 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	