

## Anti-TPD52 antibody (20-70) {FITC} (STJ503378)

STJ503378

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-TPD52 (20-70) is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.65-0.69 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:1, 000 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	7163
<b>Gene Symbol</b>	TPD52
<b>Uniprot ID</b>	TPD52_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 20-70 on human Prostate leucine zipper isoform 5 protein.
<b>Immunogen Region</b>	20-70
<b>Specificity</b>	
<b>Immunogen Sequence</b>	