

Anti-TRAF6 antibody {Biotin} (STJ503382)

STJ503382

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TRAF6 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.64-0.72 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7189
Gene Symbol	TRAF6
Uniprot ID	TRAF6_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on TRAF6 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	