

Anti-TICAM2 antibody (STJ503384)

STJ503384

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TICAM2 is suitable for use in Dot Blot, ELISA and Western Blot research applications.
Applications	DB, ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:750 DB: 1:20, 000 ELISA: 1:20, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100302736/353376
Gene Symbol	TICAM2
Uniprot ID	TCAM2_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on TRIF-related adapter molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-TICAM2 Antibody (STJ503384) with TRAM peptide. Antibody at 1:500 dilution in buffer.