

Anti-TREX1 antibody (220-245) {Biotin} (STJ503388)

STJ503388

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TREX1 (220-245) is suitable for use in ELISA, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IF, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100 IF: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11277
Gene Symbol	TREX1
Uniprot ID	TREX1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 220-245 on human TREX1 protein.
Immunogen Region	220-245
Specificity	This antibody is specific for TREX1. This antibody will not detect TREX2.
Immunogen Sequence	