

## Anti-TREX1 antibody (220-245) {FITC} (STJ503389)

STJ503389

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-TREX1 (220-245) is suitable for use in ELISA, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IHC, IF, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.55-0.75 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:10, 000 IHC: 1:50-1:100 IF: 1:100
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">11277</a>
<b>Gene Symbol</b>	<a href="#">TREX1</a>
<b>Uniprot ID</b>	<a href="#">TREX1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 220-245 on human TREX1 protein.
<b>Immunogen Region</b>	220-245
<b>Specificity</b>	This antibody is specific for TREX1. This antibody will not detect TREX2.
<b>Immunogen Sequence</b>	