

Anti-NTRK1 antibody (360-410) (STJ503405)

STJ503405

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-NTRK1 (360-410) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Chicken
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4914
Gene Symbol	NTRK1
Uniprot ID	NTRK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 360-410 on human TrkA protein.
Immunogen Region	360-410
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat brain with Anti-NTRK1 primary antibody (STJ503405) at 1:100 dilution in buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.