

Anti-TICRR antibody (STJ503417)

STJ503417

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TICRR is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	90381
Gene Symbol	TICRR
Uniprot ID	TICRR_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human Treslin protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-TICRR Antibody (STJ503417) on T47D cells. Antibody at 1:500 dilution in buffer.