

Anti-TSG101 antibody (330-390) {Biotin} (STJ503421)

STJ503421

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TSG101 (330-390) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150 ICC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7251
Gene Symbol	TSG101
Uniprot ID	TS101_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 330-390 on Tumor susceptibility gene 101 protein.
Immunogen Region	330-390
Specificity	
Immunogen Sequence	