

Anti-TTC27 antibody (100-150) {Biotin} (STJ503434)

STJ503434

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TTC27 (100-150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Zebra Fish

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74196
Gene Symbol	Ttc27
Uniprot ID	TTC27_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 100-150 on Tetratricopeptide repeat protein 27 Danio rerio protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	