

Anti-UBR1 antibody (1260-1500) {Biotin} (STJ503453)

STJ503453

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-UBR1 (1260-1500) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.70 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	197131
Gene Symbol	UBR1
Uniprot ID	UBR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1260-1500 on human UBR1 protein.
Immunogen Region	1260-1500
Specificity	
Immunogen Sequence	