

Anti-UBR1 antibody (1260-1500) {FITC} (STJ503454)

STJ503454

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-UBR1 (1260-1500) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IHC, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.68-0.70 µg/µl |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150 IHC: 1:50-1:150 |
| Formulation | |
| Isotype | IgG |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 197131 |
| Gene Symbol | UBR1 |
| Uniprot ID | UBR1_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 1260-1500 on human UBR1 protein. |
| Immunogen Region | 1260-1500 |
| Specificity | |
| Immunogen Sequence | |