

## Anti-UCN2 antibody (70-115) {Biotin} (STJ503456)

STJ503456

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-UCN2 (70-115) is suitable for use in ELISA, Dot Blot and Western Blot research applications. |
| <b>Applications</b>      | ELISA, SB, WB  |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.65-0.71 µg/µl   |
| <b>Conjugation</b>         | Biotin  |
| <b>Purification</b>        | Affinity Purified   |
| <b>Dilution Range</b>      | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000                  |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | <a href="#">90226</a>   |
| <b>Gene Symbol</b>        | <a href="#">UCN2</a>  |
| <b>Uniprot ID</b>         | <a href="#">UCN2_HUMAN</a>  |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 70-115 of human Urocortin 2 protein. |
| <b>Immunogen Region</b>   | 70-115  |
| <b>Specificity</b>        |   |
| <b>Immunogen Sequence</b> |   |