

Anti-UCN2 antibody (70-115) {FITC} (STJ503457)

STJ503457

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-UCN2 (70-115) is suitable for use in ELISA, Dot Blot and Western Blot research applications.
Applications	ELISA, SB, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.71 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	90226
Gene Symbol	UCN2
Uniprot ID	UCN2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 70-115 of human Urocortin 2 protein.
Immunogen Region	70-115
Specificity	
Immunogen Sequence	