

Anti-ULK1 antibody (1000-1050) (STJ503461)

STJ503461

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ULK1 (1000-1050) is suitable for use in ELISA, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IF, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Monkey, Mouse, Panda

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:250 ELISA: 1:10,000 IP: 1:200 IF: 1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8408
Gene Symbol	ULK1
Uniprot ID	ULK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1000-1050 on mouse ULK1 protein.
Immunogen Region	1000-1050
Specificity	
Immunogen Sequence	



Western blot with Anti-ULK1 Antibody (STJ503461) on ULK1 peptide. Antibody at 1:250 dilution in buffer.