

Anti-ULK3 antibody (452-472) (STJ503467)

STJ503467

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ULK3 (452-472) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	691171
Gene Symbol	Ulk3
Uniprot ID	ULK3_RAT
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 452-472 of Mouse ULK3, coupled to a carrier protein.
Immunogen Region	452-472
Specificity	
Immunogen Sequence	



Western blot with Anti-ULK3 Antibody (STJ503467) with ULK3 peptide. Antibody at 1:500 dilution in buffer.