

## Anti-USH2A antibody (255-305) (STJ503474)

STJ503474

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-USH2A (255-305) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IF, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.62 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution</b>	WB: 1:500
<b>Range</b>	DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">7399</a>
<b>Gene Symbol</b>	<a href="#">USH2A</a>
<b>Uniprot ID</b>	<a href="#">USH2A_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 255-305 on human Usherin isoform A protein.
<b>Immunogen Region</b>	255-305
<b>Specificity</b>	This antibody will detect both Usherin isoform A and Usherin isoform B.
<b>Immunogen Sequence</b>	