

Anti-VGluT1 antibody (1-50) {FITC} (STJ503487)

STJ503487

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VGluT1 (1-50) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	72961
Gene Symbol	Slc17a7
Uniprot ID	VGLU1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1-50 on mouse Vesicular Glutamate transporter 1 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	