

Anti-VKORC1 antibody (10-60) {FITC} (STJ503490)

STJ503490

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VKORC1 (10-60) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65-0.70 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:50, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	309004
Gene Symbol	Vkorc1
Uniprot ID	VKOR1_RAT
Immunogen	Synthetic peptide taken within amino acid region 10-60 on human VKORC1 protein.
Immunogen Region	10-60
Specificity	
Immunogen Sequence	