

Anti-VPRBP antibody (1-27) {FITC} (STJ503493)

STJ503493

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VPRBP (1-27) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	HIV

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 1-27 on human VPR protein.
Immunogen Region	1-27
Specificity	
Immunogen Sequence	