

Anti-SNF8 antibody (200-250) {Biotin} (STJ503495)

STJ503495

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SNF8 (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11267
Gene Symbol	SNF8
Uniprot ID	SNF8_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Vacuolar-sorting Protein SNF8 protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	