

Anti-NAMPT antibody (400-490) (STJ503510)

STJ503510

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NAMPT (400-490) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Pig, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.70 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:4, 000 ELISA: 1:4, 000 IP: 1:150 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10135
Gene Symbol	NAMPT
Uniprot ID	NAMPT_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 400-490 on human Visfatin protein.
Immunogen Region	400-490
Specificity	
Immunogen Sequence	



Western blot with Anti-NAMPT Antibody (STJ503510) on 1) Rat cerebral cortex, 2) Rat hippocampus, 3) Rat heart. Antibody at 1:500 dilution in buffer.