

Anti-WDR5 antibody (157-174) {Biotin} (STJ503519)

STJ503519

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WDR5 (157-174) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.67 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11091
Gene Symbol	WDR5
Uniprot ID	WDR5_HUMAN
Immunogen	Synthetic peptide: DVKTGKCLKTLPAHSDPV corresponding to amino acids 157-174 of Human WDR5 protein.
Immunogen Region	157-174
Specificity	
Immunogen Sequence	