

Anti-WEE2 antibody (60-110) (STJ503524)

STJ503524

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WEE2 (60-110) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	381759
Gene Symbol	Wee2
Uniprot ID	WEE2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 60-110 on mouse Wee1-like protein kinase isoform 2.
Immunogen Region	60-110
Specificity	
Immunogen Sequence	