

Anti-MTAP antibody (233-283) {KLH} (STJ503536)

STJ503536

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MTAP (233-283) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	KLH
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4507
Gene Symbol	MTAP
Uniprot ID	MTAP_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 233-283 on human methylthioadenosine phosphorylase protein.
Immunogen Region	233-283
Specificity	
Immunogen Sequence	