

Anti-WWOX 1 + 4 antibody (STJ503537)

STJ503537

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-WWOX 1 + 4 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.62 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000-1:100, 000 ELISA: 1:10, 000-1:100, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51741
Gene Symbol	WWOX
Uniprot ID	WWOX_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on WWOX 1 and WWOX4 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-WWOX 1 + 4 Antibody (STJ503537) with WWOX1 peptide. Antibody at 1:500 dilution in buffer.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081