

## Anti-MIEN1 antibody (40-90) {FITC} (STJ503559)

STJ503559

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-MIEN1 (40-90) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.65-0.83 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:1,000 DB: 1:2,000 ELISA: 1:2,000 IP: 1:200 IHC: 1:50-1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">84299</a>
<b>Gene Symbol</b>	<a href="#">MIEN1</a>
<b>Uniprot ID</b>	<a href="#">MIEN1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 40-90 on human XTP4 protein.
<b>Immunogen Region</b>	40-90
<b>Specificity</b>	
<b>Immunogen Sequence</b>	