

Anti-YES1 antibody (450-480) {Biotin} (STJ503565)

STJ503565

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-YES1 (450-480) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IHC, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Monkey, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.55-0.75 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100 |
| Formulation | |
| Isotype | IgG |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 7525 |
| Gene Symbol | YES1 |
| Uniprot ID | YES_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 450-480 on human YES1 protein. Host Rabbit Species Cross Reactivity Human, Monkey, Mouse, Rat |
| Immunogen Region | 450-480 |
| Specificity | |
| Immunogen Sequence | |