

Anti-YES1 antibody (450-480) {FITC} (STJ503566)

STJ503566

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-YES1 (450-480) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7525
Gene Symbol	YES1
Uniprot ID	YES_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-480 on human YES1 protein. Verified Applications ELISA, IHC, IP, WB Host Rabbit Species Cross Reactivity Human, Monkey, Mouse, Rat
Immunogen Region	450-480
Specificity	
Immunogen Sequence	