

Anti-ZFP36L1 antibody (200-250) {FITC} (STJ503575)

STJ503575

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZFP36L1 (200-250) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	677
Gene Symbol	ZFP36L1
Uniprot ID	TISB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Zinc finger protein 36, C3H1 type-like 1.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	