

## Anti-Zika Virus NS5 antibody (3230-3280) {Biotin} (STJ503592)

STJ503592

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Zika Virus NS5 (3230-3280) is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Zika

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.55-0.75 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: ELISA: 1:10, 000
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 3230-3280 on Zika virus NS5 protein.
<b>Immunogen Region</b>	3230-3280
<b>Specificity</b>	
<b>Immunogen Sequence</b>	