

## Anti-ZMPSTE24 antibody (450-480) (STJ503597)

STJ503597

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-ZMPSTE24 (450-480) is suitable for use in ELISA and Western Blot research applications. |
| <b>Applications</b>      | ELISA, WB   |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Rat, Mouse  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.55 µg/µl  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | Affinity Purified   |
| <b>Dilution Range</b>      | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000                  |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | <a href="#">230709</a>  |
| <b>Gene Symbol</b>        | <a href="#">Zmpste24</a>  |
| <b>Uniprot ID</b>         | <a href="#">FACE1_MOUSE</a>   |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 450-480 on mouse ZMPSTE24 protein. |
| <b>Immunogen Region</b>   | 450-480   |
| <b>Specificity</b>        |   |
| <b>Immunogen Sequence</b> |   |