

Anti-ZMPSTE24 antibody (450-480) {Biotin} (STJ503601)

STJ503601

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZMPSTE24 (450-480) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10269
Gene Symbol	ZMPSTE24
Uniprot ID	FACE1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-480 on mouse ZMPSTE24 protein.
Immunogen Region	450-480
Specificity	
Immunogen Sequence	