

Anti-NOS2 antibody (1064-1114) {FITC} (STJ503608)

STJ503608

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NOS2 (1064-1114) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4843
Gene Symbol	NOS2
Uniprot ID	NOS2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1064-1114 on human inducible nitric oxide synthase protein.
Immunogen Region	1064-1114
Specificity	This antibody will detect both iNOS and iNOS2 protein.
Immunogen Sequence	