

## 5HT1A protein (STJ503611)

STJ503611

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Proteins   |
| <b>Short Description</b> | 5HT1A-protein if for use in research applications. |
| <b>Applications</b>      | Immunodepletion, Immunocompetition                 |
| <b>Host/Source</b>       |  |
| <b>Reactivity</b>        |  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      |   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | 24473  |
| <b>Gene Symbol</b>        | Htr1a  |
| <b>Uniprot ID</b>         | 5HT1A_RAT  |
| <b>Immunogen</b>          | Synthetic peptide corresponding to amino acids within region 200-250 on rat 5-hydroxytryptamine receptor 1A protein. |
| <b>Immunogen Region</b>   | 200-250  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |