

5HT1B protein (STJ503614)
STJ503614**GENERAL INFORMATION**

Product Type	Proteins
Short Description	5HT1B-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	15551
Gene Symbol	Htr1b
Uniprot ID	5HT1B_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	