

TNFAIP3 protein (STJ503615)

STJ503615

GENERAL INFORMATION

Product Type	Proteins
Short Description	TNFAIP3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7128
Gene Symbol	TNFAIP3
Uniprot ID	TNAP3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 690-790 on human Tumor necrosis factor alpha-induced protein 3.
Immunogen Region	690-790
Specificity	
Immunogen Sequence	