

CCDC98 protein (STJ503617)

STJ503617

GENERAL INFORMATION

Product Type	Proteins
Short Description	CCDC98-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84142
Gene Symbol	ABRAXAS1
Uniprot ID	ABRX1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-400 on human BRCA1-A complex subunit Abraxas protein.
Immunogen Region	350-400
Specificity	
Immunogen Sequence	