

Atypical Chemokine Receptor D6 protein (STJ503622)

STJ503622

GENERAL INFORMATION

Product Type	Proteins
Short Description	Atypical Chemokine Receptor D6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1238
Gene Symbol	ACKR2
Uniprot ID	ACKR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-300 on ACRD6 protein.
Immunogen Region	200-300
Specificity	
Immunogen Sequence	